

# Prefense Monitoring FAQs

## Prefense Monitoring

### Supports Vital Signs and Waveforms for:

- ECG
- NIBP
- SpO2
- Respiratory Rate

### Monitor displays numeric vital signs and supports review of ECG

### Supports lethal arrhythmias such as:

- Asystole
- V-Fib
- V-Tach

### Alarm Events:

- Up to 9 minutes of waveform data surrounding each event will be available
- Waves and events will be stored for 72 hours

### Dedicated monitor tech 24/7:

- Monitors and reviews patient
- Retrieves tracings of changes
- Notifies primary nurse

